قسم تقنيات الاشعة المرحلة الاولى المحاضرة الخامسة فسلجه عملي



* Red blood cell count BBCS

Assist. Lecturer Doa'a Adil Almusawi

Assist. Lecturer Douaa S. Altaee



*A red blood cell count

*A red blood cell count/ is a blood test that tells you how many red blood cells (RBCs) you have.

***Other name /erythrocyte count.

***Red blood cells contain a substance called <u>haemoglobin</u>, which transports oxygen around the body.

***The amount of oxygen that's delivered to your body's tissues depends on the number of red blood cells you have

*normal range for an RBC count?

*** The normal RBC range for <u>men</u> is (<u>4.7 to 6.1</u>)million cells per microliter (mcL).

*** The normal RBC range for <u>women</u> (is 4.2 to 5.4) million/ microliter mcL.

*** The normal RBC range for <u>children</u> is (4.0 to 5.5) million/ microliter mcL.



*1. he test is almost always a part of a <u>complete bloc</u> <u>count (CBC)</u> test

*2. RBC count can be used to help diagnose bloodrelated conditions, such as <u>iron deficiency</u> <u>anaemia</u> (where there are less red blood cells than normal).

*Factors affecting R.B.C number

* Physiological factors

*1. Age, gender, activity, nutrition, pregnancy and lactation.

- * 2. In very high areas, the number of R.B.C increases due to the lack of oxygen
- * Pathological factors
- *1. Anemia.
- *2.. Leukemia
- *3. Erythropoietin hormone disorder
- *4. Polycythaemia red blood

*Polycythaemia

*Polycythaemia/ is an increase in the number of red blood cells as a result of infection of the bone marrow with cancer, which stimulates it to produce red blood cells, and its symptoms include redness of the face

POLYCYTHEMIA VERA

* INCREASED BLOOD CELL LEVELS * OVERPRODUCTION BY BONE MARROW





- *1. Haemocytometer
- *2. Red blood cells pipette
- *3. Isotonic diluting fluid(Hayem's Fluid)
- *4. Lancet:
- *5. blood
- *6. cotton
- *7.alcohol



* Hayem's Fluid)/isotonic to the Red blood cells and does not cause any damage to it.

*which preserve and fix the Red blood cells.

***The composition of Hayem's diluting Fluid:

| COMPONENTS | QUANTITY |
|-------------------|------------|
| Mercuric Chloride | 0.25 grams |
| Sodium sulfate | 2.5 grams |
| Sodium chloride | 0.5 grams |
| Distilled water | 100 ml |

*Purpose of Hayem's Fluid)

- *1. The benefit of using this solution is that it is a neutral solution that dilutes the blood and prevents its dissolution and prevents the phenomenon of sedimentation
- *. 2.Sodium sulfate/ has an effect on preventing coagulation
- *. 3.Mercury chloride/ <u>is considered as an Acts</u> as antiseptic sterilizer


pipete of red blood cell



Hayem's Fluid



* haemocytometer chamber



*procedure

*1. Take 3.98 ml of RBC diluting fluid in a Clean, Dry Test tube.

- *2. add 20 μl of Blood Specimen to the tube containing diluting fluid.
- *3. Mix well for few minutes and ready your Hemocytometer / Neubauer's Chamber
- *4. put cover slide on Neubauer's Chamber
- *5.put sample of mixing fluid on Neubauer's Chamber
- *6. under microscope count red blood cell (Calculate the five squares of the center square for counting red blood cells)
- *7. write reported



Sum of 5 squares x 10000

=Nx10000

