Practical Pharmacognosy

3^{rd.} Stage

1st semester

Dr. Zahraa Shubber

Lab.3



Saponin Glycosides

- This group of glycoside is widely distributed in higher plants.
- Saponin glycosides form **colloidal solution** in water that foam upon **shaking**,
- This is due to a decrease in the surface tension action done by saponin glycosides, as a result of the hydrophobic/ hydrophilic characteristics of the saponin, and due to this property the saponins are used in the manufacturing of beer, and soap.




- Saponins have a bitter, acrid taste, and drugs containing them are usually sternutatory and otherwise irritating the mucus membrane.
- They destroy red blood corpuscles by hemolysis and are toxic especially to cold blooded animals therefore many saponins are used as fish poisons.
- The more poisonous saponin is often called *sapotoxin*, many are toxic to insects and mollusks, and some are used to control schistosomiasis snails.

Saponin upon hydrolysis yield an aglycone known as *sapogenin*, which are crystallized upon acetylation, therefore this process is used for purification.

According to the structure of the aglycone, two kinds of saponin are recognized:

1. Pentacyclictriterpenoid saponins (acidic, and the C-atom is C30)

2. Steroidal saponins (neutral C- atom is C27).



Isolation & Identification of the Saponin Glycosides:

Procedure:

Method of extraction: Decoction.

Plant used: Saponaria officinalis family Caryophyllaceae.

Part used: Dry root.



Saponaria officinalis

Add 0.1 gm of Saponaria root in coarse powder to 20 ml distilled water in a beaker and boil gently for 2-3 minutes.

Filter hot and allow cooling

Dilute 5ml of the filtrate with water and shake vigorously.

add to the filtrate 5ml of dilute H₂SO₄ acid and boil gently for 3-5 min₅.

Filter


Supernatant is glycone

Precipitate is aglycone

^{*}The aglycones are obtained by acid hydrolysis and are insoluble in water but are soluble in 90% alcohol.




Procedure









The Chemical Tests

THE FEHLING TEST:

Aim: Identity test (specific) for Saponin glycosides.

Equipment & Reagents:

- -Test tube.
- -NaOH.
- -Fehling Reagent.
- -litmus paper.

Procedure:


- -Make the filtrate alkaline with NaOH.
- -Add 1 ml of Fehling Reagent to 3 ml of the solution.
- -heat for 10mins on boiling water bath.

Results:

Brick- red precipitation in the solution

The result of Fehling's test





Thank you

