

# Practical Pharmacognosy

3<sup>rd</sup> Stage

1<sup>st</sup> semester

**Dr. Zahraa Shubber**

**Lab.7**





# The Chemical Tests

## A. General Reaction:

### **1-Cyanidin reaction with magnesium powder.**

Flavanones and dihydroflavanols at the presence of HCl produce bright **red colour**; isoflavonoids and flavanes develop **yellow**, sometimes red colour, flavonoles do intensive **red colour**.





**2-Reaction of flavonoids with two o-oxygroups** in B cycle with lead acetate causes precipitation.

Flavones produce intensive **yellow**, aurones - **red**, anthocyanes **red or blue** precipitates.



**3-Wilson's reaction:** 5-oxyflavones and 5-oxyflavonoles with Wilson's reagent (boric and citric acids in anhydrous acetone) develop brightly **yellow** colour with **yellowish-green** fluorescence.





**4-alkaline solution:** In alkaline solutions flavanones produce **uncoloured or yellow** precipitates, that for some time become **brightly yellow or yellow** (isomerization to formation of chalcones); chalcones and aurones develop **red or purple** (it's their specific reaction), flavones and flavonoles produce **yellow coloured** precipitates.



**5-Mineral acids:** Flavones and flavonoles with mineral acids form oxonic (flavic) salts of brightly yellow or red colour, chalcones and aurones produce intensive colour of raspberry or red.



**6-Other methods:** Other methods of identification include chromatography, colorimetric or spectrophotometric analysis after reaction with aluminium chloride.

**spectrophotometric**







## **Identification of Flavonoids By Chromatography:**

### **1-By the use of Paper chromatography (P.C):**

- The stationary phase = Filter paper ( Whatman no.1).
- The mobile phase = n-BuOH:HOAc:H<sub>2</sub>O (4:1:5).
- The standard compound = Rutin.
- The spray reagent = 5% alcoholic KOH.
- Mechanism of separation = Partition.
- Developing = Ascending.




## 2-By the use of Thin layer chromatography (T.L.C):


- The stationary phase = Silica gel G.
- The mobile phase = Ethyl acetate - formic acid - glacial acetic acid - water (100:11:11:26).
- The standard compound = Rutin.
- The spray reagent = flavonoids spot on TLC plates produce a yellow-brown spots when reacted with Iodine vapor.
- Mechanism of separation = Adsorption.
- Developing = Ascending.
- Detection: Flavonoids may appear as dark spots on a green background fluoresce when observed in UV light at 254 nm. UV-plates containing fluorescent indicator (such as silica gel F254).

## Thin-layer chromatography

Step 1



Step 2



Step 3

