



Experiment 6: Melting Temperature Analysis

Department of Biochemistry

By :Abbas Hamza Khudair

Academic Year 2025-2026

Temperature Effect on Nucleic Acid Stability (Melting Temperature Analysis)

Objective

To evaluate the thermal stability of DNA duplexes and determine their melting temperature (T_m) by monitoring absorbance changes with increasing temperature.

Introduction

Thermal denaturation of DNA occurs when increasing temperature disrupts hydrogen bonding and base stacking interactions, converting double-stranded DNA (dsDNA) into single-stranded DNA (ssDNA). This transition increases absorbance at 260 nm, known as the hyperchromic effect. The melting temperature (T_m)—the point at which half of the duplex becomes denatured—provides valuable information about nucleotide sequence, GC content, length, and ionic strength. T_m analysis is central to PCR optimization, primer design, molecular hybridization, and studies of DNA–ligand interactions.

Principle

As dsDNA is heated, nucleobases become unstacked, increasing absorbance at 260 nm. Monitoring this absorbance across a temperature gradient produces a sigmoidal melting curve. T_m is derived from the curve's midpoint or from the peak of the first derivative (dA/dT). GC-rich sequences typically exhibit higher T_m values due to stronger base stacking and three hydrogen bonds between G and C.



Materials and Equipment

| Material / Equipment | Purpose |
|---|---|
| Synthetic complementary oligonucleotides | Preparation of defined duplexes |
| Tris-HCl buffer | Maintains stable pH |
| NaCl or MgCl ₂ | Modulates ionic strength and duplex stability |
| Nuclease-free water | Buffer and dilution preparation |
| Thermal cycler or heating block | Controlled temperature changes |
| UV–Vis spectrophotometer with temperature control | Monitoring absorbance vs. temperature |
| Quartz cuvettes | Heat-resistant optical measurement |
| Pipettes and filtered tips | Accurate preparation of samples |
| Microcentrifuge tubes | Duplex annealing and storage |
| Data acquisition software | Analysis of melting curves |

Procedure

1. Prepare equimolar complementary oligonucleotides (1–5 μM) in Tris buffer.
2. Denature at 95°C for 5 minutes, then allow gradual cooling to promote annealing.
3. Transfer the annealed duplex solution into a quartz cuvette.
4. Place the cuvette in the temperature-controlled UV–Vis spectrophotometer.
5. Program a heating ramp (e.g., 10°C \rightarrow 95°C at 0.5–1°C/min).
6. Record A₂₆₀ continuously or at defined temperature intervals.
7. Plot absorbance vs. temperature to produce the melting curve.
8. Determine T_m from the midpoint or derivative plot.
9. Compare T_m values under different ionic-strength conditions or sequence designs.



References

1. Bloomfield VA, Crothers DM, Tinoco I. *Nucleic Acids: Structures, Properties, and Functions*. University Science Books.
2. Watson JD et al. *Molecular Biology of the Gene*.
3. Nelson DL, Cox MM. *Lehninger Principles of Biochemistry*.
4. Cantor CR, Schimmel PR. *Biophysical Chemistry*. W.H. Freeman.
5. SantaLucia J. "DNA nearest-neighbor thermodynamics." *PNAS*.

Short Questions (No Answers)

1. How does increasing ionic strength influence DNA melting temperature?
 2. Which molecular forces are disrupted during DNA denaturation?
 3. Why do GC-rich sequences exhibit higher T_m values?
-