



Ministry of Higher education and Scientific Research

Al- Mustaqbal University

College of Science

M.S.C. SAJA JAWAD ABAID

SAJA.JAWAD.ABAID@UOMUS.EDU.IQ

LECTURE (5)

(QUALITATIVE ANALYSIS OF CARBOHYDRATE – MALTOS)

Practical Lecture on the Qualitative Analysis of Carbohydrates - Maltose

Qualitative analysis of carbohydrates involves various tests that help identify different types of carbohydrates, including disaccharides like maltose. Maltose is a disaccharide composed of two glucose molecules linked by a glycosidic bond

:Here are the steps for the qualitative analysis of maltose

1- Fehling's Test .

:Materials Needed

(Fehling's solution A (Copper(II) sulfate •

(Fehling's solution B (Alkaline tartrate solution •

:Procedure

.Add a few drops of the maltose solution into a test tube •

.Add equal volumes of Fehling's solution A and B •

.Heat the mixture gently in a boiling water bath •

If maltose is present, the solution will change from blue to red or orange due to the formation of copper(I) oxide (Cu_2O) .

: Interpretation

Maltose is a reducing sugar, which means it contains an aldehyde group capable of reducing copper(II) ions to copper(I) oxide.

2-Molisch's Test

: Materials Needed

(Molisch's reagent (α -naphthol in ethanol) •

(Concentrated sulfuric acid (H_2SO_4) •

: Procedure

Add a few drops of Molisch's reagent to the maltose solution in a test tube

Carefully add concentrated sulfuric acid along the side of the test tube to form a layer at the bottom

If a purple or violet ring forms at the interface between the acid and the solution, it indicates the presence of carbohydrates

: Interpretation

Molisch's test reacts with all carbohydrates to produce a purple or violet ring due to the formation of an aromatic compound, which confirms the presence of carbohydrates in the sample

3 - Benedict's Test .

: Materials Needed

**Benedict's reagent (a blue solution containing copper(II) •
(sulfate**

(Heat source (e.g., a boiling water bath •

: Procedure

.Add a few drops of maltose solution to a test tube •

.Add Benedict's reagent to the tube and mix •

.Heat the mixture in a boiling water bath •

**A positive result will cause the solution to change color •
from blue to green, yellow, orange, or red, depending on
.the amount of reducing sugar present**

: Interpretation

Benedict's test detects reducing sugars like maltose, which can reduce copper(II) ions to copper(I), forming a colored precipitate

4-Iodine Test (for Starch Presence)

Although this test is generally used for starch detection, it can be used to differentiate between starch and maltose

:Materials Needed

(Iodine solution (KI in water •

:Procedure

Add a few drops of iodine solution to the maltose • solution

If no blue color develops, this indicates that starch is • absent, as iodine typically forms a blue complex with starch

:Interpretation

Maltose does not react with iodine to form a blue color, differentiating it from starch

•

THANK YOU AND GOOD LUCK