




# **Blood Smear Staining**

## **4<sup>th</sup> Practical**

### **1<sup>st</sup> Term**



**Prepared and Presented by:**

**Lecturer Dr/ Ayad AbdelSalam**  
**Assist. Lecturer Dr/ Ghadeer Talib**

**Teaching of Physiology**  
**College of Technology & Health Sciences**  
**Radiology Techniques Department**

- B. Dry microscopic slide on filter paper
- C. Immerse slide in Field's stain B (Eosin) for 5 seconds
- D. Immediately wash with water
- E. Immerse slide in Field's stain A (Methylene blue) for 10 sec
- F. Immediately wash with water
- G. Dry thin films



A. Smear loopful of microbes onto slide



B. Air-dry



C. Drip methanol onto specimen to fix



D. Flood slide with stain



E. Rinse with water  
Blot dry



F. Examine with  
 $\times 100$  objective  
(oil immersion)



